



HOME THERAPY RED LIGHT SLIPPER

User Manual

Relax Muscles + Relieve Cramps

Accelerate Healing + Relieve Pain

Enhance Recovery + Circulation

before using this product read this user manual carefully and keep for reference

PACKAGE CONTENTS



1. 1x Red & Infrared LED Light Therapy Slipper (includes 66 red lights & 66 infrared lights)
2. Remote (with 20-minutes auto shut-off timer)
3. 1x AC Adapter(AC 100-240V) & Power Cord
4. 1x DC Cord

Model:	RLT-F
Material:	75% styrene butadiene rubber 25% velvet
Effective area:	12.5" × 5.3"
Power supply:	AC 100-240V or DC 12V DGYAO Portable power Bank
Frequency:	50/60 Hz
Wavelength:	660nm+880nm
Temperature :	Max 104 °F
N.W:	268g
G.W:	803g

INTRODUCTION

Thank you for purchasing your DGYAO red & infrared light therapy!

Each slipper has built-in infrared lights (880nm wavelength)+red lights (660nm wavelength). Studies have shown that a combination of infrared +red light therapy increases circulation, accelerates nerve regeneration, reduces pain, relaxes muscles, and more.

We recommend using your slippers at least once a day - in 20 minute sessions - to enjoy the full benefits of this powerful light therapy.



GETTING STARTED



Follow the below instructions to get started and enjoy the benefits of your new device:

1. Connect the slipper to either of the two DC ports on the remote control unit* Note: Remote has two DC ports. The 2nd port is a spare port for connecting a second slipper (sold separately).
2. Connect the AC Adapter to the remote
3. Connect power cord to the AC Adapter
4. Plug power cord into outlet (use USA to UK/EU adapter in the UK/Europe)
5. Press the ON/OFF switch on the remote to turn the lights on* Note: Lights will auto-shut off after 20 minutes. You may turn them back on with the ON/OFF switch. You may use the slippers multiple times a day, but do NOT exceed 40 minutes of use in one session.
6. For maximum benefit, we recommend using your light therapy slipper with bare feet (nosocks)

RED LIGHT THERAPY WARRANTY

Limited One-Year Warranty

DGYAO carries a limited one-year warranty for the original purchaser against any technical defects. Warranty coverage applies to the original owner and to the original product only and is not transferable.

DGYAO promises that the device has one year warranty for any defects in material and workmanship from the date of purchase when the device is used and maintained according to the requirements outlined in the User Manual.

If your device stops functioning or working as intended during the one-year warranty period, please call DGYAO Customer Care at Skype +86-139-2922-1675 and one of our customer service team members will work with you to understand the issue and give you instructions on how to get your DGYAO LIGHT THERAPY repaired or replaced, free of charge.

This warranty does not cover normal wear of the device and does not apply to any device that has been tampered with or damaged due to misuse, abuse, or negligent handling.

Specific Warranty Policy

- 1-Year defects liability period.
- 60 days return and replacement warranty 18-month free replacement for charger, controller and connectors. (if they are broken). Return postage beared by customers.

4.

 **HOME THERAPY RENEW**
Advanced DGYAO Support

- Return requests as “ I don’t like it”, “I don’t want it”, “the color doesn’t t look good” are not accepted.
- The situation of less than three light broken will not covered in the warranty policy because it doesn’t affect usage effect.
- Keep the original outer package for successful return.

Frequently Asked Questions

Q1: What should therapy feel like?

A1: You should feel mild warm-up on your foot and muscle relax.

Q2: How often should I use the wrap?

A2: You may use it for 1-3 times per day, 20 minutes each time. Take a rest and drink some water to relax yourself before next following use.

Q3: Does it have both red and infrared lights?

A3: Yes. 66 near infrared lights + 66 red lights.

Q4: How many lights are there on the pad?

A4: 132 lights.

Q5: What is the demension of the system?

A5: 12.5” ×5.3” Therapy area.

Q6: How long is the connector strips?

A6: Total 15 ft (4.56 meters).

Q7: Does it need to be plugged into the socket?

A7: Yes.

Q8: How long will it take to experience relief?

A8: We recommend wearing it for 2-3 therapy sessions per day for the first month. Please be patient because it may take several weeks to experience optimal relief.

Q9: Can I use this device in Europe?

A9: This device is equipped with a US plug, you need to buy a US to EU adaptor to use it.

Q10: Can I use it while eating ?

A10: Yes.

Q11: Are there side effects?

A11: No. You may experience the skin irratation if you’ re allergic to the fibers.

Q12: Can people with heart disease or heart pacemaker use it?

A12: No.

5.

 **HOME THERAPY RENEW**
Advanced DGYAO Support

Q13: What situations should I not use the therapy?

A13: When you suffer high fever, tumour, open tuberculosis, heart disease, acute soft tissue injury.

Q14: Can multiple people share the device?

A14: Yes. The device can be used by different people and your pets too if you like.

Q15: At what voltage can this device be used?

A15: The device can be used in AC 110-240V or DC 12V DGYAO Portable power Bank

Q16: What’s the wavelenth of the lights?

A16: 660nm red lights and 880nm infrared lights.

Q17: Can I put the device into water for washing?

A17: No.

Q18: Can this be used with 240 volts?

A18: Yes, It’ s fine with 100-240V.

Customer Review

★★★★★ Easy to use!

Reviewed in the United States on December 26, 2019

Color: Two Pads Set [Verified Purchase](#)

The design easy to use. I like the 20 minutes auto timer that’s built into the power button. The product is easy to set up & use. I purchased the item to help relieve my foot pain. So far, I'm thoroughly satisfied.

4 people found this helpful

[Helpful](#) [Comment](#) [Report abuse](#)

★★★★★ Amazing pain relief!

Reviewed in the United States on January 23, 2020

Color: Two Pads Set [Verified Purchase](#)

I haven't used LED therapy before and this has been amazing. Before I started using my pain was so bad it was hard to put weight in my foot. I saw results in two days. Easy to use, non-invasive and amazing results! Really helps with inflammation.

3 people found this helpful

[Helpful](#) [Comment](#) [Report abuse](#)

6.

 **HOME THERAPY RENEW**
Advanced DGYAO Support

★★★★★ Help for foot neuropathy

Reviewed in the United States on may 5, 2020

Color: One Pads Set [Verified Purchase](#)

I purchased this item after watching a webinar about neuropathy that used the red LED-infrared light therapy for treatment of neuropathy. I have been using this devise now for a couple of weeks and the neuropathy in both of my feet appears to be greatly reduced. I seems that I can start to "feel" both feet again and excited that continuing to use this devise will help even more. I was only able to purchase the single shoe/boot and hoping to get an additional one sometime in the future. I do recommend this product.

2 people found this helpful

[Helpful](#) [1 Comment](#) [Report abuse](#)

Contact Us

- Consumers and clinicians call: +86-139-2922-1675 (Skype)- for all other questions or customer care needs, please email: dgyao@dgyao.net
- Customer Care Hours: 6am to 11am, 6pm-12pm Eastern Standard Time Monday – Sunday Except major holidays.
- Manufacturers’ warranties may not apply in all cases, Please review the warranty carefully, and contact the manufacturer if you have any questions.

DONG GUAN TOP YAO INDUSTRY CO.,LTD

- ▶ Email: dgyao@dgyao.net
- ▶ Facebook: @dgyaocn
- ▶ Skype: +86 139 2922 1675
- ▶ Amazon Store: <https://www.amazon.com/dgyao>
- ▶ Official Website: <https://www.dgyao.net>



<http://www.dgyao.net>



FACEBOOK



TWITTER



YOUTUBE



LINKEDIN

7.

 **HOME THERAPY RENEW**
Advanced DGYAO Support